



## SAV-FITC (Streptavidin FITC Conjugates)

Catalog Number: FF1001

### Quantity

Flogen<sup>®</sup> SAV-FITC is provided at 1 mg/ml in sodium phosphate buffer, pH 7.2, and 0.09% of Sodium Azide.

### Requirements

Flogen<sup>®</sup> SAV-FITC can be excited at 488 nm light and leads to a fluorescence emission maximum of ~520 nm.

### Application

Flogen<sup>®</sup> SAV-FITC is a useful second-step reagent for staining with biotinylated antibodies as an avidin/biotin labeling system for Flow Cytometry, FIA and Immuno-histochemistry.

### Direction

Mixed well before use. For maximal recovery of contents, please quick-spin vial before opening.

Flogen<sup>®</sup> SAV-FITC was used as staining reagent. Since variety of applications, we recommend an appropriated dilution to be determined for further individual use.

### Storage

The Flogen<sup>®</sup> SAV-FITC should be stored in the dark at 4°C as an undiluted liquid. Do not freeze. Dilute only prior to immediate use. For research use only.